

ITGA2 / CD49b Rabbit anti-Human Polyclonal (aa1151-1180) (FITC) Antibody - LS-C260463 - LSBio

CatalogID:	LS-C260463
Target:	integrin, alpha 2 (CD49B, alpha 2 subunit of VLA-2 receptor) (ITGA2)
Synonyms:	ITGA2 Antibody, Alpha2 integrin Antibody, BR Antibody, CD49b antigen Antibody, CD49B Antibody, Integrin alpha2 Antibody, Integrin Alpha 2 Antibody, Integrin alpha-2 Antibody, HPA-5 Antibody, VLA-2 Antibody, VLAA2 Antibody, VLA-2 subunit alpha Antibody, Collagen receptor Antibody, GPIa Antibody, Platelet antigen Br Antibody, Platelet glycoprotein GPIa Antibody, VLA2 receptor, alpha-2 subunit Antibody
Family / Subfamily:	Integrin / not assigned-Integrin
Host	ITGA2 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	ITGA2 / CD49b antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	ITGA2 / CD49b antibody was raised against kLH-conjugated synthetic peptide (aa1151-1180) from the C-terminus of human ITGA2.
Specificity:	Human ITGA2 / CD49b
Epitope:	aa1151-1180
Reactivity:	Human, Mouse
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/29/2014

© 2014 LifeSpan BioSciences