

ITGB2 / MAC-1 / CD18 Rabbit anti-Human Polyclonal (FITC) Antibody - LS-C260513 - LSBio

CatalogID:	LS-C260513
Target:	integrin, beta 2 (complement component 3 receptor 3 and 4 subunit) (ITGB2)
Synonyms:	ITGB2 Antibody, CD18 antigen Antibody, CD18 Antibody, Integrin beta chain, beta 2 Antibody, LCAMB Antibody, MAC-1 Antibody, MF17 Antibody, MF17 Antibody, Integrin beta-2 Antibody, LAD Antibody, LFA-1 Antibody
Family / Subfamily:	Integrin / not assigned-Integrin
Host	ITGB2 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	ITGB2 / MAC-1 / CD18 antibody was raised against Human
Antigen Type:	Recombinant protein
Immunogen:	ITGB2 / MAC-1 / CD18 antibody was raised against recombinant human ITGB2.
Specificity:	Human ITGB2 / MAC-1 / CD18
Reactivity:	Human
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/29/2014

© 2014 LifeSpan BioSciences