

ITPKB Rabbit anti-Human Polyclonal (aa31-61) (PE) Antibody - LS-C260580 - LSBio

CatalogID:	LS-C260580
Target:	inositol-trisphosphate 3-kinase B (ITPKB)
Synonyms:	ITPKB Antibody, IP3 3-kinase B Antibody, IP3K-B Antibody, IP3-3KB Antibody, IP3KB Antibody, InsP 3-kinase B Antibody, IP3K Antibody, IP3K B Antibody, PIG37 Antibody
Family / Subfamily:	Non-protein Kinase / ITPK
Host	ITPKB antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	ITPKB antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	ITPKB antibody was raised against kLH-conjugated synthetic peptide (aa31-61) from the N-terminal region of human ITPKB.
Specificity:	Human ITPKB.
Epitope:	aa31-61
Reactivity:	Human
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, Western blot, Flow Cytometry, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/29/2014

© 2014 LifeSpan BioSciences