

KDR / VEGFR2 Rabbit anti-Human Polyclonal (Tyr1175) (FITC) Antibody - LS-C261048 - LSBio	
CatalogID:	LS-C261048
Target:	kinase insert domain receptor (a type III receptor tyrosine kinase) (KDR)
Synonyms:	KDR Antibody, CD309 Antibody, Fetal liver kinase 1 Antibody, Fetal liver kinase-1 Antibody, FLK-1 Antibody, FLK1 Antibody, Kinase insert domain receptor Antibody, Quek1 Antibody, VEGFR Antibody, VEGFR2 Antibody, Soluble VEGFR2 Antibody, VEGF receptor 2 Antibody, VEGFR-2 Antibody, CD309 antigen Antibody, VEGFR2 Antibody
Family / Subfamily:	Protein Kinase / VEGF Receptor
Host	KDR antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	KDR / VEGFR2 antibody was raised against Human
Antigen Type:	Synthetic phospho-peptide
Immunogen:	KDR / VEGFR2 antibody was raised against kLH-conjugated synthetic phosphopeptide (aa1160-1190) around Tyr1175 of human KDR.
Specificity:	Human KDR / VEGFR2
Epitope:	Tyr1175
Reactivity:	Human
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	ELISA, Dot Blot (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/29/2014 © 2014 LifeSpan BioSciences	