

KHNYN Rabbit anti-Human Polyclonal (aa586-615) (FITC) Antibody - LS-C261080 - LSBio	
CatalogID:	LS-C261080
Target:	KH and NYN domain containing (KHNYN)
Synonyms:	KHNYN Antibody, KH and NYN domain containing Antibody, KIAA0323 Antibody, Protein KHNYN Antibody
Host	KHNYN antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	KHNYN antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	KHNYN antibody was raised against kLH-conjugated synthetic peptide (aa586-615) from the C-terminus of human KHNYN.
Specificity:	Human KHNYN
Epitope:	aa586-615
Reactivity:	Human
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only Not for resale without prior written consent from LifeSpan BioSciences, Inc. Created on 9/29/2014 © 2014 LifeSpan BioSciences	