

| C-Kit / CD117 Rabbit anti-Human Polyclonal (aa688-718) (FITC) Antibody - LS-C261194 - LSBio |   |
|---|---|
| <b>CatalogID:</b>   | LS-C261194  |
| <b>Target:</b>  | v-kit Hardy-Zuckerman 4 feline sarcoma viral oncogene homolog   |
| <b>Synonyms:</b>  | KIT Antibody, C-Kit Antibody, CD117 Antibody, PBT Antibody, Piebald trait Antibody, Piebald trait protein Antibody, Proto-oncogene c-Kit Antibody, SCFR Antibody, Soluble KIT variant 1 Antibody, Tyrosine-protein kinase Kit Antibody, p145 c-kit Antibody, CD117 antigen Antibody |
| <b>Family / Subfamily:</b>  | Protein Kinase / PDGF Receptor  |
| <b>Host</b>   | KIT antibody was produced in Rabbit   |
| <b>Clonality:</b>   | Polyclonal  |
| <b>Isotype:</b>   | IgG   |
| <b>Conjugations:</b>  | Fluorescein (FITC)  |
| <b>Immunogen Species:</b>   | C-Kit / CD117 antibody was raised against Human   |
| <b>Antigen Type:</b>  | Synthetic peptide   |
| <b>Immunogen:</b>   | C-Kit / CD117 antibody was raised against kLH-conjugated synthetic peptide (aa688-718) from the Y703 region of human KIT.   |
| <b>Specificity:</b>   | Human C-Kit / CD117   |
| <b>Epitope:</b>   | aa688-718   |
| <b>Reactivity:</b>  | Human, Mouse  |
| <b>Purification:</b>  | Affinity purified   |
| <b>Presentation:</b>  | PBS, pH 7.2, 0.09% sodium azide.  |
| <b>Recommended Storage:</b>   | May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.   |
| <b>Usage Summary:</b>   | The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.  |
| <b>Uses:</b>  | Western blot, Flow Cytometry, ELISA (Optimal dilution to be determined by the researcher)   |
| <b>Size:</b>  | 200 µl  |
| <b>Requested From:</b>  | Japan   |
| <b>Laboratory Reagent For In Vitro Research Use Only</b>                                    |   |
| <b>Not for resale without prior written consent from LifeSpan BioSciences, Inc.</b>         |   |
| Created on 9/29/2014<br>© 2014 LifeSpan BioSciences   |   |