

| C-Kit / CD117 Rabbit anti-Human Polyclonal (aa807-836) (FITC) Antibody - LS-C261199 - LSBio | |
|---|---|
| CatalogID: | LS-C261199 |
| Target: | v-kit Hardy-Zuckerman 4 feline sarcoma viral oncogene homolog |
| Synonyms: | KIT Antibody, C-Kit Antibody, CD117 Antibody, PBT Antibody, Piebald trait Antibody, Piebald trait protein Antibody, Proto-oncogene c-Kit Antibody, SCFR Antibody, Soluble KIT variant 1 Antibody, Tyrosine-protein kinase Kit Antibody, p145 c-kit Antibody, CD117 antigen Antibody |
| Family / Subfamily: | Protein Kinase / PDGF Receptor |
| Host | KIT antibody was produced in Rabbit |
| Clonality: | Polyclonal |
| Isotype: | IgG |
| Conjugations: | Fluorescein (FITC) |
| Immunogen Species: | C-Kit / CD117 antibody was raised against Human |
| Antigen Type: | Synthetic peptide |
| Immunogen: | C-Kit / CD117 antibody was raised against kLH-conjugated synthetic peptide (aa807-836) from the C-terminus of human KIT. |
| Specificity: | Human C-Kit / CD117 |
| Epitope: | aa807-836 |
| Reactivity: | Human |
| Purification: | Affinity purified |
| Presentation: | PBS, pH 7.2, 0.09% sodium azide. |
| Recommended Storage: | May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months. |
| Usage Summary: | The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested. |
| Uses: | Western blot, ELISA (Optimal dilution to be determined by the researcher) |
| Size: | 200 µl |
| Requested From: | Japan |
| Laboratory Reagent For In Vitro Research Use Only | |
| Not for resale without prior written consent from LifeSpan BioSciences, Inc. | |
| Created on 9/29/2014 © 2014 LifeSpan BioSciences | |