

C-Kit / CD117 Rabbit anti-Human Polyclonal (pSer821) (PE) Antibody - LS-C261298 - LSBio

CatalogID:	LS-C261298
Target:	v-kit Hardy-Zuckerman 4 feline sarcoma viral oncogene homolog
Synonyms:	KIT Antibody, C-Kit Antibody, CD117 Antibody, PBT Antibody, Piebald trait Antibody, Piebald trait protein Antibody, Proto-oncogene c-Kit Antibody, SCFR Antibody, Soluble KIT variant 1 Antibody, Tyrosine-protein kinase Kit Antibody, p145 c-kit Antibody, CD117 antigen Antibody
Family / Subfamily:	Protein Kinase / PDGF Receptor
Host	KIT antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	C-Kit / CD117 antibody was raised against Human
Antigen Type:	Synthetic phospho-peptide
Immunogen:	C-Kit / CD117 antibody was raised against kLH-conjugated synthetic phosphopeptide around Ser821 of human KIT.
Specificity:	Human C-Kit / CD117
Epitope:	pSer821
Reactivity:	Human
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	ELISA, Dot Blot (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/29/2014

© 2014 LifeSpan BioSciences