

C-Kit / CD117 Rabbit anti-Human Polyclonal (pThr718) (FITC) Antibody - LS-C261332 - LSBio

| | |
|-----------------------------|---|
| CatalogID: | LS-C261332 |
| Target: | v-kit Hardy-Zuckerman 4 feline sarcoma viral oncogene homolog |
| Synonyms: | KIT Antibody, C-Kit Antibody, CD117 Antibody, PBT Antibody, Piebald trait Antibody, Piebald trait protein Antibody, Proto-oncogene c-Kit Antibody, SCFR Antibody, Soluble KIT variant 1 Antibody, Tyrosine-protein kinase Kit Antibody, p145 c-kit Antibody, CD117 antigen Antibody |
| Family / Subfamily: | Protein Kinase / PDGF Receptor |
| Host | KIT antibody was produced in Rabbit |
| Clonality: | Polyclonal |
| Isotype: | IgG |
| Conjugations: | Fluorescein (FITC) |
| Immunogen Species: | C-Kit / CD117 antibody was raised against Human |
| Antigen Type: | Synthetic phospho-peptide |
| Immunogen: | C-Kit / CD117 antibody was raised against kLH-conjugated synthetic phosphopeptide around Thr718 of human KIT. |
| Specificity: | Human C-Kit / CD117 |
| Epitope: | pThr718 |
| Reactivity: | Human |
| Purification: | Affinity purified |
| Presentation: | PBS, pH 7.2, 0.09% sodium azide. |
| Recommended Storage: | May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months. |
| Usage Summary: | The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested. |
| Uses: | ELISA, Dot Blot (Optimal dilution to be determined by the researcher) |
| Size: | 200 µl |
| Requested From: | Japan |

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/29/2014

© 2014 LifeSpan BioSciences