

C-Kit / CD117 Rabbit anti-Mouse Polyclonal (pTyr719) (FITC) Antibody - LS-C261368 - LSBio	
CatalogID:	LS-C261368
Target:	v-kit Hardy-Zuckerman 4 feline sarcoma viral oncogene homolog
Synonyms:	KIT Antibody, C-Kit Antibody, CD117 Antibody, PBT Antibody, Piebald trait Antibody, Piebald trait protein Antibody, Proto-oncogene c-Kit Antibody, SCFR Antibody, Soluble KIT variant 1 Antibody, Tyrosine-protein kinase Kit Antibody, p145 c-kit Antibody, CD117 antigen Antibody
Family / Subfamily:	Protein Kinase / PDGF Receptor
Host	KIT antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	C-Kit / CD117 antibody was raised against Mouse
Antigen Type:	Synthetic phospho-peptide
Immunogen:	C-Kit / CD117 antibody was raised against kLH-conjugated synthetic phosphopeptide around Tyr719 of mouse KIT.
Specificity:	Mouse C-Kit / CD117
Epitope:	pTyr719
Reactivity:	Mouse
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	ELISA, Dot Blot (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/29/2014	
© 2014 LifeSpan BioSciences	