

C-Kit / CD117 Rabbit anti-Mouse Polyclonal (pTyr719) (PE) Antibody - LS-C261370 - LSBio	
<b>CatalogID:</b>	LS-C261370
<b>Target:</b>	v-kit Hardy-Zuckerman 4 feline sarcoma viral oncogene homolog
<b>Synonyms:</b>	KIT Antibody, C-Kit Antibody, CD117 Antibody, PBT Antibody, Piebald trait Antibody, Piebald trait protein Antibody, Proto-oncogene c-Kit Antibody, SCFR Antibody, Soluble KIT variant 1 Antibody, Tyrosine-protein kinase Kit Antibody, p145 c-kit Antibody, CD117 antigen Antibody
<b>Family / Subfamily:</b>	Protein Kinase / PDGF Receptor
<b>Host</b>	KIT antibody was produced in Rabbit
<b>Clonality:</b>	Polyclonal
<b>Isotype:</b>	IgG
<b>Conjugations:</b>	Phycoerythrin (PE)
<b>Immunogen Species:</b>	C-Kit / CD117 antibody was raised against Mouse
<b>Antigen Type:</b>	Synthetic phospho-peptide
<b>Immunogen:</b>	C-Kit / CD117 antibody was raised against kLH-conjugated synthetic phosphopeptide around Tyr719 of mouse KIT.
<b>Specificity:</b>	Mouse C-Kit / CD117
<b>Epitope:</b>	pTyr719
<b>Reactivity:</b>	Mouse
<b>Purification:</b>	Affinity purified
<b>Presentation:</b>	PBS, pH 7.2, 0.09% sodium azide.
<b>Recommended Storage:</b>	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
<b>Usage Summary:</b>	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
<b>Uses:</b>	ELISA, Dot Blot (Optimal dilution to be determined by the researcher)
<b>Size:</b>	200 µl
<b>Requested From:</b>	Japan
<b>Laboratory Reagent For In Vitro Research Use Only</b>	
<b>Not for resale without prior written consent from LifeSpan BioSciences, Inc.</b>	
Created on 9/29/2014 © 2014 LifeSpan BioSciences	