

KLF4 Rabbit anti-Human Polyclonal (Ser245) (AP) Antibody - LS-C261425 - LSBio

CatalogID:	LS-C261425
Target:	Kruppel-like factor 4 (gut) (KLF4)
Synonyms:	KLF4 Antibody, EZF Antibody, GKLf Antibody, Kruppel-like factor 4 (gut) Antibody, Hezf Antibody, Krueppel-like factor 4 Antibody
Host	KLF4 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Alkaline Phosphatase (AP)
Immunogen Species:	KLF4 antibody was raised against Human
Antigen Type:	Synthetic phospho-peptide
Immunogen:	KLF4 antibody was raised against kLH-conjugated synthetic phosphopeptide (aa132-161) around Ser245 of human KLF4.
Specificity:	Human KLF4.
Epitope:	Ser245
Reactivity:	Human
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	ICC, ELISA, Dot Blot (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/29/2014

© 2014 LifeSpan BioSciences