

LAMP2 / CD107b Rabbit anti-Human Polyclonal (FITC) Antibody - LS-C261679 - LSBio

| | |
|-----------------------------|---|
| CatalogID: | LS-C261679 |
| Target: | lysosomal-associated membrane protein 2 (LAMP2) |
| Synonyms: | LAMP2 Antibody, CD107b Antibody, CD107b antigen Antibody, LAMPB Antibody, LGP110 Antibody, LAMP-2 Antibody |
| Family / Subfamily: | Cancer |
| Host | LAMP2 antibody was produced in Rabbit |
| Clonality: | Polyclonal |
| Isotype: | IgG |
| Conjugations: | Fluorescein (FITC) |
| Immunogen Species: | LAMP2 / CD107b antibody was raised against Human |
| Antigen Type: | Recombinant protein |
| Immunogen: | LAMP2 / CD107b antibody was raised against recombinant Human LAMP2. |
| Specificity: | Human LAMP2 / CD107b. |
| Reactivity: | Human |
| Purification: | Affinity purified |
| Presentation: | PBS, pH 7.2, 0.09% sodium azide. |
| Recommended Storage: | May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months. |
| Usage Summary: | The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested. |
| Uses: | Western blot, ELISA (Optimal dilution to be determined by the researcher) |
| Size: | 200 µl |
| Requested From: | Japan |

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/29/2014

© 2014 LifeSpan BioSciences