

LAMP2 / CD107b Rabbit anti-Human Polyclonal (aa203-233) (FITC) Antibody - LS-C261691 - LSBio

CatalogID:	LS-C261691
Target:	lysosomal-associated membrane protein 2 (LAMP2)
Synonyms:	LAMP2 Antibody, CD107b Antibody, CD107b antigen Antibody, LAMPB Antibody, LGP110 Antibody, LAMP-2 Antibody
Family / Subfamily:	Cancer
Host	LAMP2 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	LAMP2 / CD107b antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	LAMP2 / CD107b antibody was raised against kLH-conjugated synthetic peptide (aa203-233) from internal region of human LAMP2.
Specificity:	Human LAMP2 / CD107b.
Epitope:	aa203-233
Reactivity:	Human
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/29/2014

© 2014 LifeSpan BioSciences