

LIMK1 / LIMK Rabbit anti-Human Polyclonal (Thr508) (PE) Antibody - LS-C261922 - LSBio	
CatalogID:	LS-C261922
Target:	LIM domain kinase 1 (LIMK1)
Synonyms:	LIMK1 Antibody, LIM domain kinase 1 Antibody, Lim kinase 1 Antibody, LIMK-1 Antibody, LIM kinase Antibody, LIMK Antibody
Family / Subfamily:	Protein Kinase / LIMK
Host	LIMK1 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	LIMK1 / LIMK antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	LIMK1 / LIMK antibody was raised against kLH-conjugated synthetic peptide selected surrounding Thr508 of human LIMK1.
Specificity:	Human LIMK1 / LIMK
Epitope:	Thr508
Reactivity:	Human
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	ELISA, Dot Blot (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/29/2014 © 2014 LifeSpan BioSciences	