

LRTM2 Rabbit anti-Human Polyclonal (aa346-370) (PE) Antibody - LS-C262312 - LSBio

CatalogID:	LS-C262312
Target:	leucine-rich repeats and transmembrane domains 2 (LRTM2)
Synonyms:	LRTM2 Antibody
Host	LRTM2 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	LRTM2 antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	LRTM2 antibody was raised against kLH-conjugated synthetic peptide (aa346-370) from the C-terminus of human LRTM2.
Specificity:	Human LRTM2
Epitope:	aa346-370
Reactivity:	Human
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, ICC, Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/29/2014

© 2014 LifeSpan BioSciences