

**LYN Rabbit anti-Mouse Polyclonal (aa1-30) (FITC) Antibody - LS-C262388 - LSBio**

<b>CatalogID:</b>	LS-C262388
<b>Target:</b>	v-yes-1 Yamaguchi sarcoma viral related oncogene homolog (LYN)
<b>Synonyms:</b>	LYN Antibody, Lyn protein tyrosine kinase Antibody, Lyn ptk Antibody, Lyn tyrosine kinase Antibody, p56Lyn Antibody, p53Lyn Antibody, Tyrosine-protein kinase Lyn Antibody, JTK8 Antibody, Lyn kinase Antibody
<b>Family / Subfamily:</b>	Protein Kinase / Src
<b>Host</b>	LYN antibody was produced in Rabbit
<b>Clonality:</b>	Polyclonal
<b>Isotype:</b>	IgG
<b>Conjugations:</b>	Fluorescein (FITC)
<b>Immunogen Species:</b>	LYN antibody was raised against Mouse
<b>Antigen Type:</b>	Synthetic peptide
<b>Immunogen:</b>	LYN antibody was raised against kLH-conjugated synthetic peptide (aa1-30) from the N-terminal region of Mouse Lyn.
<b>Specificity:</b>	Mouse LYN
<b>Epitope:</b>	aa1-30
<b>Reactivity:</b>	Mouse
<b>Purification:</b>	Affinity purified
<b>Presentation:</b>	PBS, pH 7.2, 0.09% sodium azide.
<b>Recommended Storage:</b>	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
<b>Usage Summary:</b>	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
<b>Uses:</b>	IHC, Western blot, ELISA (Optimal dilution to be determined by the researcher)
<b>Size:</b>	200 µl
<b>Requested From:</b>	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/29/2014

© 2014 LifeSpan BioSciences