

LYN Rabbit anti-Mouse Polyclonal (aa1-30) (PE) Antibody - LS-C262392 - LSBio

CatalogID:	LS-C262392
Target:	v-yes-1 Yamaguchi sarcoma viral related oncogene homolog (LYN)
Synonyms:	LYN Antibody, Lyn protein tyrosine kinase Antibody, Lyn ptk Antibody, Lyn tyrosine kinase Antibody, p56Lyn Antibody, p53Lyn Antibody, Tyrosine-protein kinase Lyn Antibody, JTK8 Antibody, Lyn kinase Antibody
Family / Subfamily:	Protein Kinase / Src
Host	LYN antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	LYN antibody was raised against Mouse
Antigen Type:	Synthetic peptide
Immunogen:	LYN antibody was raised against kLH-conjugated synthetic peptide (aa1-30) from the N-terminal region of Mouse Lyn.
Specificity:	Mouse LYN
Epitope:	aa1-30
Reactivity:	Mouse
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/29/2014

© 2014 LifeSpan BioSciences