

MAP1LC3A / LC3A Rabbit anti-Human Polyclonal (Thr6) (FITC) Antibody - LS-C262507 - LSBio	
CatalogID:	LS-C262507
Target:	microtubule-associated protein 1 light chain 3 alpha (MAP1LC3A)
Synonyms:	MAP1LC3A Antibody, ATG8E Antibody, LC3 Antibody, LC3A Antibody, MAP1BLC3 Antibody, MAP1A/1B light chain 3 A Antibody, MAP1A/MAP1B LC3 A Antibody, MAP1A/MAP1B light chain 3 A Antibody, MAP1ALC3 Antibody
Host	MAP1LC3A antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	MAP1LC3A / LC3A antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	MAP1LC3A / LC3A antibody was raised against kLH-conjugated synthetic peptide (aa1-30) around Thr6 of human MAP1LC3B.
Specificity:	Human MAP1LC3A / LC3A.
Epitope:	Thr6
Reactivity:	Human, Mouse
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/29/2014	
© 2014 LifeSpan BioSciences	