

**MAP1LC3A / LC3A Rabbit anti-Human Polyclonal (aa81-111) (FITC) Antibody - LS-C262516 - LSBio**

<b>CatalogID:</b>	LS-C262516
<b>Target:</b>	microtubule-associated protein 1 light chain 3 alpha (MAP1LC3A)
<b>Synonyms:</b>	MAP1LC3A Antibody, ATG8E Antibody, LC3 Antibody, LC3A Antibody, MAP1BLC3 Antibody, MAP1A/1B light chain 3 A Antibody, MAP1A/MAP1B LC3 A Antibody, MAP1A/MAP1B light chain 3 A Antibody, MAP1ALC3 Antibody
<b>Host</b>	MAP1LC3A antibody was produced in Rabbit
<b>Clonality:</b>	Polyclonal
<b>Isotype:</b>	IgG
<b>Conjugations:</b>	Fluorescein (FITC)
<b>Immunogen Species:</b>	MAP1LC3A / LC3A antibody was raised against Human
<b>Antigen Type:</b>	Synthetic peptide
<b>Immunogen:</b>	MAP1LC3A / LC3A antibody was raised against kLH-conjugated synthetic peptide (aa81-111) from internal region of human PMAT-Y198.
<b>Specificity:</b>	Human MAP1LC3A / LC3A.
<b>Epitope:</b>	aa81-111
<b>Reactivity:</b>	Human, Mouse
<b>Purification:</b>	Affinity purified
<b>Presentation:</b>	PBS, pH 7.2, 0.09% sodium azide.
<b>Recommended Storage:</b>	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
<b>Usage Summary:</b>	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
<b>Uses:</b>	IHC, Western blot, Flow Cytometry, ELISA (Optimal dilution to be determined by the researcher)
<b>Size:</b>	200 µl
<b>Requested From:</b>	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/29/2014

© 2014 LifeSpan BioSciences