

MECP2 Rabbit anti-Human Polyclonal (pSer80) (PE) Antibody - LS-C262927 - LSBio

CatalogID:	LS-C262927
Target:	methyl CpG binding protein 2 (Rett syndrome) (MECP2)
Synonyms:	MECP2 Antibody, AUTSX3 Antibody, Methyl-CpG-binding protein 2 Antibody, MRX16 Antibody, MRX79 Antibody, MRXSL Antibody, RTT Antibody, RTS Antibody, MeCp-2 protein Antibody, RS Antibody, MRXS13 Antibody, PPMX Antibody
Host	MECP2 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	MECP2 antibody was raised against Human
Antigen Type:	Synthetic phospho-peptide
Immunogen:	MECP2 antibody was raised against kLH-conjugated synthetic phosphopeptide around Ser80 of human MECP2.
Specificity:	Human MECP2.
Epitope:	pSer80
Reactivity:	Human
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	ELISA, Dot Blot (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/29/2014

© 2014 LifeSpan BioSciences