

METTL2A Rabbit anti-Human Polyclonal (aa336-367) (FITC) Antibody - LS-C263070 - LSBio

CatalogID:	LS-C263070
Target:	Methyltransferase Like 2A (METTL2A)
Synonyms:	METTL2A Antibody, Methyltransferase like 2A Antibody, METTL2 Antibody
Host	METTL2A antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	METTL2A antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	METTL2A antibody was raised against kLH-conjugated synthetic peptide (aa336-367) from the C-terminus of human METTL2.
Specificity:	Human METTL2A.
Epitope:	aa336-367
Reactivity:	Human
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, ICC, Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/29/2014

© 2014 LifeSpan BioSciences