

MMP9 / Gelatinase B Rabbit anti-Human Polyclonal (aa548-578) (FITC) Antibody - LS-C263322 - LSBio

CatalogID:	LS-C263322
Target:	matrix metalloproteinase 9 (gelatinase B, 92kDa gelatinase, 92kDa type IV collagenase) (MMP9)
Synonyms:	MMP9 Antibody, 92kd gelatinase Antibody, 92 kd type iv collagenase Antibody, 92 kDa gelatinase Antibody, 92 kDa type IV collagenase Antibody, Matrix metalloproteinase-9 Antibody, MMP-9 Antibody, Type V collagenase Antibody, CLG4B Antibody, Gelatinase B Antibody, GELB Antibody, Macrophage gelatinase Antibody, MANDP2 Antibody
Family / Subfamily:	Protease / Metalloproteinase M10A
Host	MMP9 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	MMP9 / Gelatinase B antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	MMP9 / Gelatinase B antibody was raised against kLH-conjugated synthetic peptide (aa548-578) from the C-terminus of human MMP9.
Specificity:	Human MMP9 / MMP-9
Epitope:	aa548-578
Reactivity:	Human
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/29/2014

© 2014 LifeSpan BioSciences