

MUSK Rabbit anti-Human Polyclonal (aa42-72) (FITC) Antibody - LS-C263757 - LSBio

CatalogID:	LS-C263757
Target:	muscle, skeletal, receptor tyrosine kinase (MUSK)
Synonyms:	MUSK Antibody, Nsk1 Antibody, Nsk2 Antibody, MDK4 Antibody, Muscle-specific MUSK receptor Antibody
Family / Subfamily:	Protein Kinase / MUSK
Host	MUSK antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	MUSK antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	MUSK antibody was raised against kLH-conjugated synthetic peptide (aa42-72) from the N-terminal region of human MUSK.
Specificity:	Human MUSK
Epitope:	aa42-72
Reactivity:	Human
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, Western blot, Flow Cytometry, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/29/2014

© 2014 LifeSpan BioSciences