

MYBPHL Rabbit anti-Human Polyclonal (aa87-116) (PE) Antibody - LS-C263813 - LSBio

CatalogID:	LS-C263813
Target:	myosin binding protein H-like (MYBPHL)
Synonyms:	MYBPHL Antibody, Myosin binding protein H-like Antibody, Myosin-binding protein H-like Antibody
Host	MYBPHL antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	MYBPHL antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	MYBPHL antibody was raised against kLH-conjugated synthetic peptide (aa87-116) from the N-terminal region of human MYBPHL.
Specificity:	Human MYBPHL
Epitope:	aa87-116
Reactivity:	Human
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, ICC, Western blot, Flow Cytometry, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/29/2014

© 2014 LifeSpan BioSciences