

MYCBP Rabbit anti-Mouse Polyclonal (aa1-30) (AP) Antibody - LS-C263814 - LSBio

CatalogID:	LS-C263814
Target:	MYC binding protein (MYCBP)
Synonyms:	MYCBP Antibody, AMY-1 Antibody, Associate of Myc 1 Antibody, C-Myc-binding protein Antibody, Associate of myc-1 Antibody, C-myc binding protein Antibody, AMY1 Antibody
Host	MYCBP antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Alkaline Phosphatase (AP)
Immunogen Species:	MYCBP antibody was raised against Mouse
Antigen Type:	Synthetic peptide
Immunogen:	MYCBP antibody was raised against kLH-conjugated synthetic peptide (aa1-30) from the N-terminal region of human MYCBP.
Specificity:	Mouse MYCBP
Epitope:	aa1-30
Reactivity:	Mouse
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/29/2014

© 2014 LifeSpan BioSciences