

NANOG Rabbit anti-Human Polyclonal (Ser285) (APC) Antibody - LS-C263947 - LSBio

CatalogID:	LS-C263947
Target:	Nanog homeobox
Synonyms:	NANOG Antibody, Nanog homeobox Antibody, NANOGP8 Antibody, HNanog Antibody, Homeobox protein NANOG Antibody
Host	NANOG antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Allophycocyanin (APC)
Immunogen Species:	NANOG antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	NANOG antibody was raised against kLH-conjugated synthetic peptide (aa272-302) around Ser285 of human NANOG.
Specificity:	Human NANOG.
Epitope:	Ser285
Reactivity:	Human
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	ICC, Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/29/2014

© 2014 LifeSpan BioSciences