

NPHS1 / Nephtrin Rabbit anti-Human Polyclonal (Tyr1193) (FITC) Antibody - LS-C264211 - LSBio	
CatalogID:	LS-C264211
Target:	nephrosis 1, congenital, Finnish type (nephrin) (NPHS1)
Synonyms:	NPHS1 Antibody, Nephtrin Antibody, NPHN Antibody, CNF Antibody
Family / Subfamily:	Immunoglobulin / not assigned-Immunoglobulin
Host	NPHS1 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	NPHS1 / Nephtrin antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	NPHS1 / Nephtrin antibody was raised against kLH-conjugated synthetic peptide selected around Tyr1193 of human Nephtrin.
Specificity:	Human Nephtrin
Epitope:	Tyr1193
Reactivity:	Human
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	ELISA, Dot Blot (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/29/2014	
© 2014 LifeSpan BioSciences	