

NLK Rabbit anti-Human Polyclonal (aa268-298) (FITC) Antibody - LS-C264482 - LSBio

| | |
|-----------------------------|---------------------------------------------------------------------------------------------------------------------------------------------|
| CatalogID: | LS-C264482 |
| Target: | nemo-like kinase (NLK) |
| Synonyms: | NLK Antibody, LAK1 Antibody, Nemo-like kinase Antibody, Nmo kinase Antibody, Nemo like kinase Antibody, Protein LAK1 Antibody |
| Family / Subfamily: | Protein Kinase / MAPK |
| Host | NLK antibody was produced in Rabbit |
| Clonality: | Polyclonal |
| Isotype: | IgG |
| Conjugations: | Fluorescein (FITC) |
| Immunogen Species: | NLK antibody was raised against Human |
| Antigen Type: | Synthetic peptide |
| Immunogen: | NLK antibody was raised against kLH-conjugated synthetic peptide (aa268-298) from internal region of human PMAT-Y198. |
| Specificity: | Human NLK |
| Epitope: | aa268-298 |
| Reactivity: | Human |
| Purification: | Affinity purified |
| Presentation: | PBS, pH 7.2, 0.09% sodium azide. |
| Recommended Storage: | May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months. |
| Usage Summary: | The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested. |
| Uses: | ICC, Western blot, Flow Cytometry, ELISA (Optimal dilution to be determined by the researcher) |
| Size: | 200 µl |
| Requested From: | Japan |

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/29/2014

© 2014 LifeSpan BioSciences