

NPC / NPC1 Rabbit anti-Human Polyclonal (aa590-620) (FITC) Antibody - LS-C264604 - LSBio	
CatalogID:	LS-C264604
Target:	Niemann-Pick disease, type C1 (NPC1)
Synonyms:	NPC1 Antibody, NPC Antibody, Niemann-Pick C1 protein Antibody, Niemann-Pick disease, type C1 Antibody
Host	NPC1 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	NPC / NPC1 antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	NPC / NPC1 antibody was raised against KLH-conjugated synthetic peptide (aa590-620) from internal region of human NPC1.
Specificity:	Human NPC / NPC1
Epitope:	aa590-620
Reactivity:	Human
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, ICC, Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/29/2014	
© 2014 LifeSpan BioSciences	