

OAF Rabbit anti-Human Polyclonal (aa54-83) (FITC) Antibody - LS-C264905 - LSBio

CatalogID:	LS-C264905
Target:	OAF homolog (Drosophila)
Synonyms:	OAF Antibody, NS5ATP13TP2 Antibody, OAF homolog (Drosophila) Antibody, Out at first protein homolog Antibody
Host	OAF antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	OAF antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	OAF antibody was raised against kLH-conjugated synthetic peptide (aa54-83) from the N-terminal region of human OAF.
Specificity:	Human OAF
Epitope:	aa54-83
Reactivity:	Human
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/29/2014

© 2014 LifeSpan BioSciences