

IL1A / IL1 Alpha Mouse Monoclonal (FITC) Antibody - LS-C26491 - LSBio

CatalogID:	LS-C26491
Target:	interleukin 1, alpha (IL1A)
Synonyms:	IL1A Antibody, Hematopoietin-1 Antibody, IL-1 alpha Antibody, IL1 Antibody, IL1-ALPHA Antibody, IL1 Alpha Antibody, Interleukin-1 alpha Antibody, IL-1A Antibody, IL1F1 Antibody, Pro-interleukin-1-alpha Antibody, Preinterleukin 1 alpha Antibody, Interleukin 1, alpha Antibody
Family / Subfamily:	IL-1 / not assigned-IL-1
Host	IL1A antibody was produced in Mouse
Clonality:	Monoclonal
Conjugations:	Fluorescein (FITC)
Reactivity:	Human
Presentation:	PBS, 0.08% Sodium Azide
Uses:	Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	100 tst
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/29/2014

© 2014 LifeSpan BioSciences