

POU5F1 / OCT4 Rabbit anti-Human Polyclonal (aa5-34) (AP) Antibody - LS-C264938 - LSBio

<b>CatalogID:</b>	LS-C264938
<b>Target:</b>	POU class 5 homeobox 1 (POU5F1)
<b>Synonyms:</b>	POU5F1 Antibody, Oct-4 Antibody, OCT3 Antibody, OCT4 Antibody, OTF4 Antibody, POU class 5 homeobox 1 Antibody, Octamer-binding protein 3 Antibody, Octamer-binding protein 4 Antibody, OTF-3 Antibody, Oct-3 Antibody, OTF3 Antibody
<b>Host</b>	POU5F1 antibody was produced in Rabbit
<b>Clonality:</b>	Polyclonal
<b>Isotype:</b>	IgG
<b>Conjugations:</b>	Alkaline Phosphatase (AP)
<b>Immunogen Species:</b>	POU5F1 / OCT4 antibody was raised against Human
<b>Antigen Type:</b>	Synthetic peptide
<b>Immunogen:</b>	POU5F1 / OCT4 antibody was raised against kLH-conjugated synthetic peptide (aa5-34) from the N-terminal region of human 39724.
<b>Specificity:</b>	Human POU5F1 / OCT4.
<b>Epitope:</b>	aa5-34
<b>Reactivity:</b>	Human
<b>Purification:</b>	Affinity purified
<b>Presentation:</b>	PBS, pH 7.2, 0.09% sodium azide.
<b>Recommended Storage:</b>	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
<b>Usage Summary:</b>	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
<b>Uses:</b>	ICC, Western blot, ELISA (Optimal dilution to be determined by the researcher)
<b>Size:</b>	200 µl
<b>Requested From:</b>	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/29/2014

© 2014 LifeSpan BioSciences