

IL1B Mouse Monoclonal (FITC) Antibody - LS-C26494 - LSBio

| | |
|----------------------------|---|
| CatalogID: | LS-C26494 |
| Target: | interleukin 1, beta (IL1B) |
| Synonyms: | IL1B Antibody, IL-1B Antibody, IL-1 beta Antibody, IL1-BETA Antibody, Preinterleukin 1 beta Antibody, Pro-interleukin-1-beta Antibody, Catabolin Antibody, IL-1 Antibody, IL1F2 Antibody, Interleukin 1, beta Antibody, Interleukin-1 beta Antibody |
| Family / Subfamily: | IL-1 / not assigned-IL-1 |
| Host | IL1B antibody was produced in Mouse |
| Clonality: | Monoclonal |
| Isotype: | IgG1 |
| Conjugations: | Fluorescein (FITC) |
| Reactivity: | Human |
| Presentation: | PBS, 0.08% sodium azide, 0.2% carrier protein, sterile-filtered |
| Uses: | Flow Cytometry (Optimal dilution to be determined by the researcher) |
| Size: | 100 tst |
| Requested From: | Japan |

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/29/2014

© 2014 LifeSpan BioSciences