

IL1RN Mouse Monoclonal (FITC) Antibody - LS-C26500 - LSBio

CatalogID:	LS-C26500
Target:	interleukin 1 receptor antagonist (IL1RN)
Synonyms:	IL1RN Antibody, Anakinra Antibody, DIRA Antibody, IL-1ra Antibody, IL-1RN Antibody, IL1F3 Antibody, IL1RA Antibody, MVCD4 Antibody, IL1 inhibitor Antibody, IL1RN (IL1F3) Antibody, IRAP Antibody, ICIL-1RA Antibody, IL-1ra3 Antibody
Family / Subfamily:	IL-1 / not assigned-IL-1
Host	IL1RN antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG1
Conjugations:	Fluorescein (FITC)
Reactivity:	Human
Presentation:	PBS, 0.08% Sodium Azide
Uses:	Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	100 tst
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/29/2014

© 2014 LifeSpan BioSciences