

IL6 Mouse Monoclonal (FITC) Antibody - LS-C26501 - LSBio

CatalogID:	LS-C26501
Target:	interleukin 6 (interferon, beta 2) (IL6)
Synonyms:	IL6 Antibody, B-cell stimulatory factor 2 Antibody, BSF-2 Antibody, BSF2 Antibody, CTL differentiation factor Antibody, Hybridoma growth factor Antibody, IFN-beta-2 Antibody, Interleukin-6 Antibody, Interferon beta-2 Antibody, Interleukin BSF-2 Antibody, IFNB2 Antibody, B-cell differentiation factor Antibody, CDF Antibody, HSF Antibody, IL-6 Antibody
Family / Subfamily:	Hormone / not assigned-Hormone
Host	IL6 antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG1
Conjugations:	Fluorescein (FITC)
Reactivity:	Human
Presentation:	PBS, 1.7% glycerol, 0.5% gelatin, 0.05% ProClin.
Uses:	Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	100 tst
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/29/2014

© 2014 LifeSpan BioSciences