

MYC / c-Myc Mouse Monoclonal Antibody - LS-C26522 - LSBio

CatalogID:	LS-C26522
Target:	v-myc avian myelocytomatosis viral oncogene homolog
Synonyms:	MYC Antibody, BHLHe39 Antibody, C-Myc Antibody, MRTL Antibody, Myc proto-oncogene protein Antibody, Proto-oncogene c-Myc Antibody, Rats1 Antibody, Transcription factor p64 Antibody, v-myc Antibody
Host	MYC antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG1
Reactivity:	Human, Avian
Presentation:	PBS, 0.08% Sodium Azide
Uses:	IHC, Western blot, Immunoprecipitation (Optimal dilution to be determined by the researcher)
Size:	100 µg
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/29/2014

© 2014 LifeSpan BioSciences