

OVOL1 Rabbit anti-Human Polyclonal (aa1-30) (APC) Antibody - LS-C265339 - LSBio

CatalogID:	LS-C265339
Target:	ovo-like zinc finger 1 (OVOL1)
Synonyms:	OVOL1 Antibody, HOVO1 Antibody, Ovo-like 1(Drosophila) Antibody, Ovo homolog-like 1 Antibody
Host	OVOL1 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Allophycocyanin (APC)
Immunogen Species:	OVOL1 antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	OVOL1 antibody was raised against kLH-conjugated synthetic peptide (aa1-30) from the N-terminal region of human OVOL1.
Specificity:	Human OVOL1
Epitope:	aa1-30
Reactivity:	Human
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/29/2014

© 2014 LifeSpan BioSciences