

OVOL1 Rabbit anti-Human Polyclonal (aa1-30) (Azide-free) (HRP) Antibody - LS-C265342 - LSBio	
CatalogID:	LS-C265342
Target:	ovo-like zinc finger 1 (OVOL1)
Synonyms:	OVOL1 Antibody, HOVO1 Antibody, Ovo-like 1(Drosophila) Antibody, Ovo homolog-like 1 Antibody
Host	OVOL1 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Horseradish Peroxidase (HRP)
Modifications:	Azide-free
Immunogen Species:	OVOL1 antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	OVOL1 antibody was raised against kLH-conjugated synthetic peptide (aa1-30) from the N-terminal region of human OVOL1.
Specificity:	Human OVOL1
Epitope:	aa1-30
Reactivity:	Human
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, azide-free.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/29/2014	
© 2014 LifeSpan BioSciences	