

POU5F1B Rabbit anti-Human Polyclonal (aa217-247) (AP) Antibody - LS-C265419 - LSBio	
CatalogID:	LS-C265419
Target:	POU class 5 homeobox 1B (POU5F1B)
Synonyms:	POU5F1B Antibody, OCT4-PG1 Antibody, POU 5 domain protein Antibody, POU5F1P1 Antibody, POU5FLC20 Antibody, POU class 5 homeobox 1B Antibody, POU5F1P4 Antibody, POU5FLC8 Antibody, Octamer-binding protein 3-like Antibody, OTF3C Antibody, OTF3P1 Antibody
Host	POU5F1B antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Alkaline Phosphatase (AP)
Immunogen Species:	POU5F1B antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	POU5F1B antibody was raised against kLH-conjugated synthetic peptide (aa217-247) from the C-terminus of human P5F1B.
Specificity:	Human POU5F1B
Epitope:	aa217-247
Reactivity:	Human
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	Western blot, Flow Cytometry, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only Not for resale without prior written consent from LifeSpan BioSciences, Inc. Created on 9/29/2014 © 2014 LifeSpan BioSciences	