

| POU5F1B Rabbit anti-Human Polyclonal (aa217-247) (PE) Antibody - LS-C265424 - LSBio | |
|---|---|
| CatalogID: | LS-C265424 |
| Target: | POU class 5 homeobox 1B (POU5F1B) |
| Synonyms: | POU5F1B Antibody, OCT4-PG1 Antibody, POU 5 domain protein Antibody, POU5F1P1 Antibody, POU5FLC20 Antibody, POU class 5 homeobox 1B Antibody, POU5F1P4 Antibody, POU5FLC8 Antibody, Octamer-binding protein 3-like Antibody, OTF3C Antibody, OTF3P1 Antibody |
| Host | POU5F1B antibody was produced in Rabbit |
| Clonality: | Polyclonal |
| Isotype: | IgG |
| Conjugations: | Phycoerythrin (PE) |
| Immunogen Species: | POU5F1B antibody was raised against Human |
| Antigen Type: | Synthetic peptide |
| Immunogen: | POU5F1B antibody was raised against kLH-conjugated synthetic peptide (aa217-247) from the C-terminus of human P5F1B. |
| Specificity: | Human POU5F1B |
| Epitope: | aa217-247 |
| Reactivity: | Human |
| Purification: | Affinity purified |
| Presentation: | PBS, pH 7.2, 0.09% sodium azide. |
| Recommended Storage: | May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months. |
| Usage Summary: | The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested. |
| Uses: | Western blot, Flow Cytometry, ELISA (Optimal dilution to be determined by the researcher) |
| Size: | 200 µl |
| Requested From: | Japan |
| Laboratory Reagent For In Vitro Research Use Only | |
| Not for resale without prior written consent from LifeSpan BioSciences, Inc. | |
| Created on 9/29/2014 | |
| © 2014 LifeSpan BioSciences | |