

PARP1 / PARP Rabbit anti-Rat Polyclonal (pSer373) (PE) Antibody - LS-C265594 - LSBio	
<b>CatalogID:</b>	LS-C265594
<b>Target:</b>	poly (ADP-ribose) polymerase 1 (PARP1)
<b>Synonyms:</b>	PARP1 Antibody, ADPRT Antibody, ADP-ribosyltransferase NAD(+) Antibody, Adp-ribosyltransferase Antibody, ADPRT1 Antibody, ARTD1 Antibody, ADPRT 1 Antibody, Poly(ADP-ribose) synthetase Antibody, Poly(ADP-ribosyl)transferase Antibody, PARP Antibody, Poly [ADP-ribose] polymerase 1 Antibody, PADPRT-1 Antibody, PARP-1 Antibody, Poly (ADP-ribose) polymerase 1 Antibody, Poly(ADP-ribose) polymerase Antibody, Poly[ADP-ribose] synthase 1 Antibody, PPOL Antibody
<b>Host</b>	PARP1 antibody was produced in Rabbit
<b>Clonality:</b>	Polyclonal
<b>Isotype:</b>	IgG
<b>Conjugations:</b>	Phycoerythrin (PE)
<b>Immunogen Species:</b>	PARP1 / PARP antibody was raised against Rat
<b>Antigen Type:</b>	Synthetic phospho-peptide
<b>Immunogen:</b>	PARP1 / PARP antibody was raised against kLH-conjugated synthetic phosphopeptide around Ser373 of rat PARP1.
<b>Specificity:</b>	Rat PARP1 / PARP
<b>Epitope:</b>	pSer373
<b>Reactivity:</b>	Rat, Rabbit
<b>Purification:</b>	Affinity purified
<b>Presentation:</b>	PBS, pH 7.2, 0.09% sodium azide.
<b>Recommended Storage:</b>	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
<b>Usage Summary:</b>	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
<b>Uses:</b>	ELISA, Dot Blot (Optimal dilution to be determined by the researcher)
<b>Size:</b>	200 µl
<b>Requested From:</b>	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/29/2014	
© 2014 LifeSpan BioSciences	