

| PARP1 / PARP Rabbit anti-Human Polyclonal (pThr368) (PE) Antibody - LS-C265618 - LSBio | |
|--|---|
| CatalogID: | LS-C265618 |
| Target: | poly (ADP-ribose) polymerase 1 (PARP1) |
| Synonyms: | PARP1 Antibody, ADPRT Antibody, ADP-ribosyltransferase NAD(+) Antibody, Adp-ribosyltransferase Antibody, ADPRT1 Antibody, ARTD1 Antibody, ADPRT 1 Antibody, Poly(ADP-ribose) synthetase Antibody, Poly(ADP-ribosyl)transferase Antibody, PARP Antibody, Poly [ADP-ribose] polymerase 1 Antibody, PADPRT-1 Antibody, PARP-1 Antibody, Poly (ADP-ribose) polymerase 1 Antibody, Poly(ADP-ribose) polymerase Antibody, Poly[ADP-ribose] synthase 1 Antibody, PPOL Antibody |
| Host | PARP1 antibody was produced in Rabbit |
| Clonality: | Polyclonal |
| Isotype: | IgG |
| Conjugations: | Phycoerythrin (PE) |
| Immunogen Species: | PARP1 / PARP antibody was raised against Human |
| Antigen Type: | Synthetic phospho-peptide |
| Immunogen: | PARP1 / PARP antibody was raised against kLH-conjugated synthetic phosphopeptide around Thr368 of human PARP1. |
| Specificity: | Human PARP1 / PARP |
| Epitope: | pThr368 |
| Reactivity: | Human |
| Purification: | Affinity purified |
| Presentation: | PBS, pH 7.2, 0.09% sodium azide. |
| Recommended Storage: | May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months. |
| Usage Summary: | The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested. |
| Uses: | ELISA, Dot Blot (Optimal dilution to be determined by the researcher) |
| Size: | 200 µl |
| Requested From: | Japan |
| Laboratory Reagent For In Vitro Research Use Only | |
| Not for resale without prior written consent from LifeSpan BioSciences, Inc. | |
| Created on 9/29/2014 | |
| © 2014 LifeSpan BioSciences | |