

PARP1 / PARP Rabbit anti-Human Polyclonal (pThr373) (PE) Antibody - LS-C265625 - LSBio	
CatalogID:	LS-C265625
Target:	poly (ADP-ribose) polymerase 1 (PARP1)
Synonyms:	PARP1 Antibody, ADPRT Antibody, ADP-ribosyltransferase NAD(+) Antibody, Adp-ribosyltransferase Antibody, ADPRT1 Antibody, ARTD1 Antibody, ADPRT 1 Antibody, Poly(ADP-ribose) synthetase Antibody, Poly(ADP-ribosyl)transferase Antibody, PARP Antibody, Poly [ADP-ribose] polymerase 1 Antibody, PADPRT-1 Antibody, PARP-1 Antibody, Poly (ADP-ribose) polymerase 1 Antibody, Poly(ADP-ribose) polymerase Antibody, Poly[ADP-ribose] synthase 1 Antibody, PPOL Antibody
Host	PARP1 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	PARP1 / PARP antibody was raised against Human
Antigen Type:	Synthetic phospho-peptide
Immunogen:	PARP1 / PARP antibody was raised against kLH-conjugated synthetic phosphopeptide around Thr373 of human PARP1.
Specificity:	Human PARP1 / PARP
Epitope:	pThr373
Reactivity:	Human
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	ELISA, Dot Blot (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/29/2014	
© 2014 LifeSpan BioSciences	