

PDGFRA / PDGFR Alpha Rabbit anti-Mouse Polyclonal (aa677-706) (FITC) Antibody - LS-C265932 - LSBio

CatalogID:	LS-C265932
Target:	platelet-derived growth factor receptor, alpha polypeptide (PDGFRA)
Synonyms:	PDGFRA Antibody, CD140A Antibody, CD140a antigen Antibody, PDGFalphaR Antibody, PDGFR2 Antibody, PDGFRalpha Antibody, AlphaPDGFR Antibody, RHEPDGFRA Antibody, PDGF receptor alpha Antibody, PDGFR-2 Antibody, PDGFRA/BCR fusion Antibody, APDGFR Antibody, Pdgf alpha-receptor Antibody, PDGF-R-alpha Antibody, PDGFR-alpha Antibody
Family / Subfamily:	Protein Kinase / PDGF Receptor
Host	PDGFRA antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	PDGFRA / PDGFR Alpha antibody was raised against Mouse
Antigen Type:	Synthetic peptide
Immunogen:	PDGFRA / PDGFR Alpha antibody was raised against kLH-conjugated synthetic peptide (aa677-706) from internal region of Mouse Pdgfra.
Specificity:	Mouse PDGFRA / PDGFR Alpha
Epitope:	aa677-706
Reactivity:	Mouse
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/29/2014

© 2014 LifeSpan BioSciences