

PDGFRA / PDGFR Alpha Rabbit anti-Human Polyclonal (pTyr762) (PE) Antibody - LS-C265944 - LSBio	
CatalogID:	LS-C265944
Target:	platelet-derived growth factor receptor, alpha polypeptide (PDGFRA)
Synonyms:	PDGFRA Antibody, CD140A Antibody, CD140a antigen Antibody, PDGFalphaR Antibody, PDGFR2 Antibody, PDGFRalpha Antibody, AlphaPDGFR Antibody, RHEPDGFRA Antibody, PDGF receptor alpha Antibody, PDGFR-2 Antibody, PDGFRA/BCR fusion Antibody, APDGFR Antibody, Pdgf alpha-receptor Antibody, PDGF-R-alpha Antibody, PDGFR-alpha Antibody
Family / Subfamily:	Protein Kinase / PDGF Receptor
Host	PDGFRA antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	PDGFRA / PDGFR Alpha antibody was raised against Human
Antigen Type:	Synthetic phospho-peptide
Immunogen:	PDGFRA / PDGFR Alpha antibody was raised against kLH-conjugated synthetic phosphopeptide around Tyr762 of human PDGFRA.
Specificity:	Human PDGFRA / PDGFR Alpha
Epitope:	pTyr762
Reactivity:	Human
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	ELISA, Dot Blot (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/29/2014	
© 2014 LifeSpan BioSciences	