

PFDN5 / MM1 Rabbit anti-Human Polyclonal (aa78-108) (PE) Antibody - LS-C266085 - LSBio

CatalogID:	LS-C266085
Target:	prefoldin subunit 5 (PFDN5)
Synonyms:	PFDN5 Antibody, C-Myc-binding protein Mm-1 Antibody, C-myc binding protein Antibody, MM-1 Antibody, Myc modulator 1 Antibody, Myc modulator-1 Antibody, Prefoldin 5 Antibody, MM1 Antibody, PFD5 Antibody, Prefoldin subunit 5 Antibody
Host	PFDN5 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	PFDN5 / MM1 antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	PFDN5 / MM1 antibody was raised against kLH-conjugated synthetic peptide (aa78-108) from internal region of human PFDN5.
Specificity:	Human PFDN5 / MM1
Epitope:	aa78-108
Reactivity:	Human, Mouse
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/29/2014

© 2014 LifeSpan BioSciences