

PIK3C3 / VPS34 Rabbit anti-Human Polyclonal (Ser676) (FITC) Antibody - LS-C266244 - LSBio	
<b>CatalogID:</b>	LS-C266244
<b>Target:</b>	phosphatidylinositol 3-kinase, catalytic subunit type 3 (PIK3C3)
<b>Synonyms:</b>	PIK3C3 Antibody, HVps34 Antibody, Phosphatidylinositol 3-kinase Antibody, VPS34 Antibody, PI3-kinase type 3 Antibody, PI3K type 3 Antibody, PtdIns-3-kinase type 3 Antibody
<b>Family / Subfamily:</b>	Non-Protein Kinase / not assigned-Non-Protein Kinase
<b>Host</b>	PIK3C3 antibody was produced in Rabbit
<b>Clonality:</b>	Polyclonal
<b>Isotype:</b>	IgG
<b>Conjugations:</b>	Fluorescein (FITC)
<b>Immunogen Species:</b>	PIK3C3 / VPS34 antibody was raised against Human
<b>Antigen Type:</b>	Synthetic peptide
<b>Immunogen:</b>	PIK3C3 / VPS34 antibody was raised against kLH-conjugated synthetic peptide (aa661-691) around Ser676 of human PI3KC3.
<b>Specificity:</b>	Human PIK3C3 / VPS34.
<b>Epitope:</b>	Ser676
<b>Reactivity:</b>	Human
<b>Purification:</b>	Affinity purified
<b>Presentation:</b>	PBS, pH 7.2, 0.09% sodium azide.
<b>Recommended Storage:</b>	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
<b>Usage Summary:</b>	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
<b>Uses:</b>	IHC, ELISA (Optimal dilution to be determined by the researcher)
<b>Size:</b>	200 µl
<b>Requested From:</b>	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/29/2014	
© 2014 LifeSpan BioSciences	