

PIK3C3 / VPS34 Rabbit anti-Human Polyclonal (aa32-61) (FITC) Antibody - LS-C266268 - LSBio	
CatalogID:	LS-C266268
Target:	phosphatidylinositol 3-kinase, catalytic subunit type 3 (PIK3C3)
Synonyms:	PIK3C3 Antibody, HVps34 Antibody, Phosphatidylinositol 3-kinase Antibody, VPS34 Antibody, PI3-kinase type 3 Antibody, PI3K type 3 Antibody, PtdIns-3-kinase type 3 Antibody
Family / Subfamily:	Non-Protein Kinase / not assigned-Non-Protein Kinase
Host	PIK3C3 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	PIK3C3 / VPS34 antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	PIK3C3 / VPS34 antibody was raised against kLH-conjugated synthetic peptide (aa32-61) from the N-terminal region of human autophagy kinase PI3KC3.
Specificity:	Human PIK3C3 / VPS34.
Epitope:	aa32-61
Reactivity:	Human
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, ICC, Western blot, Immunoprecipitation, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/29/2014	
© 2014 LifeSpan BioSciences	