

PRAMEF6 Rabbit anti-Human Polyclonal (aa402-431) (FITC) Antibody - LS-C267048 - LSBio	
CatalogID:	LS-C267048
Target:	PRAME family member 6 (PRAMEF6)
Synonyms:	PRAMEF6 Antibody, PRAME family member 6 Antibody
Host	PRAMEF6 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	PRAMEF6 antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	PRAMEF6 antibody was raised against kLH-conjugated synthetic peptide (aa402-431) from the C-terminus of human PRAMEF6.
Specificity:	Human PRAMEF6
Epitope:	aa402-431
Reactivity:	Human
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, ICC, Western blot, Flow Cytometry, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/30/2014	
© 2014 LifeSpan BioSciences	